

Nitrosamines Analysis in Pharmaceuticals

Using Triple Quadrupole LC/MS/MS and Quadrupole Time of Flight LC/MS Consumables workflow ordering guide





Mutagenic impurities in APIs and drug products pose a significant risk to health and safety even in small quantities—and thus are a major concern for drug makers. Mutagenic impurities can damage DNA, leading to mutations and potentially cancer. Efforts to address and control the presence of trace levels of mutagenic impurities is of special concern to global regulators. As a result, the U.S. FDA and other regulatory agencies have taken steps to address the issue of mutagenic impurities in pharmaceuticals¹. Detection and quantification of these trace nitrosamines in APIs and drug products can be challenging, necessitating the use of advanced and sensitive tools to meet regulatory requirements.

The list of APIs and drug products for nitrosamine determination has expanded beyond angiotensin II receptor blocker (ARB) Sartan drugs and include metformin, an oral diabetes drug and histamine-2 receptor antagonists such as ranitidine. This is evidenced by the recent recalls of metformin by various regulatory bodies like the U.S. Food & Drug Administration (FDA), European Directorate for the Quality of Medicines (EDQM), and Health Sciences Authority (HSA) due to the presence of N-nitroso-dimethylamine (NDMA). These impurities: (N-nitrosodimethylamine (NDMA), N-nitrosodiethylamine (NDEA), N-nitrosodiisopropylamine (NDIPA), N-nitrosoethylisopropylamine (NEIPA) and N-nitrosodibutylamine (NDBA) are classified as probable human carcinogens and are believed to have been introduced into the finished products due to chemical reactions that occur during the API manufacturing process.

These impurities can be detected using either a single quadrupole GC/MS (GC/MSD), or triple quadrupole GC/MS/MS (GC/TQ)(1), a triple quadrupole LC/MS/MS (LC/TQ), or quadrupole time of flight LC/MS (LC/Q-TOF)(2-7). LC/MS/MS-based methods are generally very specific and highly sensitive. For this reason, these have served as the basis for development of methods to detect and quantify nitrosamine impurities in drug substance and drug products such as metformin, valsartan, losartan and irbesartan.









Figure 1. From left to right: 1260 Infinity II LC System, 6470B Triple Quadrupole LC/MS, 6550 iFunnel Q-TOF LC/MS, and Ultivo Triple Quadruple LC/MS.

Column Choices

The pentafluorophenyl (PFP) ligand on the **InfinityLab Poroshell 120 PFP column** used in Method 2 (Table 1) (*8*) provides an orthogonal separation mechanism with C18 chemistries. PFP phases can separate analytes based on small differences in structure, substitution, and steric access to polar moieties. The resulting selectivity for positional isomers, halogenated compounds, and polar analytes is particularly useful when analyzing complex mixtures. Since NDIPA and NDPA are positional isomers, the InfinityLab Poroshell 120 PFP is the ideal recommended column for this separation for *ease of use*.

The **InfinityLab Poroshell HPH-C18 column**, which has also been used for this analysis in Method 1 (Table 1) (8) incorporates hybrid particle technology which improves particle *ruggedness at extended pH*, enabling long lifetimes and fewer column changes. When using this column, however, a robust method development process is critical to ensure that the method is long lasting, stable and reliable. Because the retention and selectivity of ionizable compounds can change significantly with varying pH, it is becoming standard practice to employ low-, medium-, and high-pH analyses during method development. In addition to optimizing the gradient conditions, in order to achieve separation between the positional isomers NDIPA and NDPA, the instrument MS/MS parameters also need to be optimized to maximize sensitivity.

LC configuration and parameters

Table 1. UHPLC configuration and settings. For method details, see Reference 8.

Parameter	Value							
	Method 1			Method 2				
Instruments	Agilent 1 (G7120A Agilent 1 (G7167E Agilent 1 mostat (Agilent 1 (G13150	290 Infinit) 290 Infinit) 290 Infinit G7116B) 260 Infinit ;)	y II high∹ y II multi y II multi y diode a	speed pump sampler column ther- array detector	Agilent 1290 Infinity II high-speed pump (G7120A) Agilent 1290 Infinity II multisampler (G7167B) Agilent 1290 Infinity II multicolumn thermostat (G7116B) Agilent 1260 Infinity diode array detector (G1315C)			
Needle Wash	Methano	ol: water (8	0:20)		Methanol: water (80:20)			
Sample Diluent	Water: methanol (95:5)				Methanol			
Multisampler Temp.	10 °C				10 °C			
Injection Volume	20 µL				5 μL			
Analytical Column	Agilent InfinityLab Poroshell HPH-C18, 4.6 × 150 mm, 2.7 μm (p/n 693975-702(T))				Agilent InfinityLab Poroshell 120 PFP, 3.0 x 150 mm, 2.7 μm (p/n 693975-308)			
Column Temp.	40 °C				40 °C			
Mobile Phase A	0.1% formic acid in water				0.1% formic acid in water			
Mobile Phase B	0.1% formic acid in methanol				0.1% formic acid in methanol			
Flow Rate	0.5 mL/min				0.5 mL/min			
Gradient	Time (min) 0 2 7 10 11 16.5 16.6 20.0	% A 95 95 40 25 10 10 95 95	% B 5 5 60 75 90 90 5 5	Flow (mL/min) 0.5 0.5 0.5 0.5 0.5 0.5 0.5 0.5 0.5	Time (min) 0 3 14 17 19 19.1 22	% A 95 95 40 10 10 95 95	% B 5 5 60 90 90 5 5	Flow (mL/min) 0.5 0.5 0.5 0.5 0.5 0.5 0.5 0.5
Stop Time	20 minutes				22 minutes			
UV Wavelength	230 nm				230 nm			



Figure 2. Representative MRM chromatogram of all the nitrosamine impurities at 0.5 ng/ mL using Method 1 (For method details see Reference 8).



Figure 3. Representative MRM chromatogram of all eight nitrosamine impurities at 3 ng/ mL using Method 2 (For method details see Reference 8).

Easy Selection and Ordering Information

This guide provides recommendations for Agilent products used in this analysis, so you can find what you're looking for quickly. Click the MyList* links in the header below to add items to your "Favorite Products" list at the Agilent online store. Then, enter the quantities for the products you need. Your list will remain under "Favorite Products" for your use with future orders.

Click MyList to download all items in this table to your Favorite Products list

Product Description	Part Number
Nitrosamine Standards and Solvents	
Nitrosamine standards (in methylene chloride)	<u>US-113N-1</u>
InfinityLab Ultrapure LCMS MeOH	<u>5191-4497</u>
InfinityLab Ultrapure LCMS water	<u>5191-4498</u>
Formic acid, 5 mL	<u>G2453-85060</u>
Nitrosamine - Sample Preparation	
15 mL Centrifuge tubes and caps, 50/pk	<u>5610-2039</u>
Captiva Econofilter, PVDF membrane, 13 mm diameter, 0.2 µm pore size, 1000/pk	<u>5190-5261</u>
Captiva Premium Syringe Filter PES membrane, 15 mm diameter, 0.2 μm pore size, 100/pk (LCMS certified)	<u>5190-5096</u> 1
Captiva Premium Syringe Filter Nylon membrane 15mm diameter, 0.2 μm pore size, 100/pk (HPLC certified)	<u>5190-5088</u> ²
Captiva Disposable Syringe, 5 mL, 100/pk	<u>9301-6476</u>
Captiva Disposable Syringe, 10 mL, 100/pk	<u>9301-6474</u>
Nitrosamine - LC Column	
InfinityLab Poroshell 120 PFP, 3.0 x 150 mm, 2.7 μm (recommended)	<u>693975-308</u>
InfinityLab Poroshell 120 PFP guard column, 3.0 mm, 2.7 μm (recommended)	<u>823750-915</u>
InfinityLab Poroshell HPH-C18, 4.6 × 150 mm, 2.7 µm	<u>693975-702</u>
InfinityLab Poroshell HPH-C18, 4.6 mm, 2.7 µm, UHPLC guard, 3/pk	<u>820750-922</u>
Nitrosamine - LC Supplies	
InfinityLab Quick Connect assembly, 0.12 x 105 mm, for column inlet connection on UHPLC	<u>5067-5957</u>
InfinityLab Quick Connect assembly, 0.17 x 105 mm, for column inlet connection on UHPLC	<u>5067-6166</u>
InfinityLab Quick Turn fitting, for column outlet	<u>5067-5966</u>
Quick Turn capillary 0.12 x 280 mm, for connection from column to detector	<u>5500-1191</u>
Kit of stay safe waste cap GL45 with 4 ports and waste can 6 L	<u>5043-1221</u>
Charcoal filter with time strip for waste container	<u>5043-1193</u>
InfinityLab Stay Safe cap, starter kit	<u>5043-1222</u>
Stainless steel solvent inlet filter, 10 µm pore size	01018-60025
InfinityLab solvent filtration assembly includes glass funnel, 250 mL, membrane holder glass base, glass flask, 1 L, and aluminum clamp	<u>5191-6776</u> ³
Regenerated cellulose membrane 47 mm, 0.20 µm 100/pk	<u>5191-4340</u> ³











Nitrosamines-Vials and Caps	
Vial, screw top, amber, write-on spot, certified, 2 mL, 100/pk. Vial size: 12 x 32 mm (12 mm cap)	<u>5182-0716</u>
Cap, screw, green, preslit PTFE/silicone, 100/pk. Cap size: 12 mm	<u>5183-2077</u>
Vial insert, 250 µL, deactivated glass with polymer feet, 100/pk 5	<u>5181-8872</u>
Nitrosamines-MS Supplies	
APCI Needle Replacement Kit	<u>G1946-68704</u>
APCI Needle	<u>G1960-20030</u>
Capillary, Fast Switching, 0.6 mm	<u>G1960-80060</u>



* First time using "MyList"? You will be asked to enter your email address for account verification. If you have an existing Agilent account, you will be able to log in. If you don't have a registered Agilent account, you will need to <u>register for one</u>. This feature is valid only in countries that are e-commerce enabled. All items can also be ordered through your regular sales and distributor channels.

1. Having similar performance to PVDF, this membrane has low protein binding and is ideal for protein analysis. Available in a smaller pack size (100/pk) and is LC/MS certified.

2. Ideal for general filtration needs. Should not be used for protein analysis. Available in a smaller pack size (100/pk) and is HPLC certified.

3. Solvent filtration assembly and associated filter membranes are not recommended for use with InfinityLab Ultrapure LC/MS solvents.

References

- 1. Nitrosamines analysis in pharmaceuticals using single quadrupole GC/MS and triple quadrupole GC/MS/MS: Consumables workflow ordering guide, Agilent publication (5994-2979EN)
- 2. Nitrosamine Impurities Application Guide Confidently detect and quantify mutagenic impurities in APIs and Drug Products (<u>5994-2393EN</u>)
- 3. Determination of a Genotoxic NDMA Impurity Using the High-Resolution Agilent 6546 LC/Q-TOF in Ranitidine Drug Substance and Drug Products (5994-1626EN)
- 4. Simultaneous Determination of Eight Nitrosamine Impurities in Metformin Using the Agilent 6470 Triple Quadrupole LC/MS (<u>5994-2286EN</u>)
- 5. Determination of Nitrosamine Impurities Using the Ultivo Triple Quadrupole LC/MS (5994-1383EN)
- 6. Determination of NDMA Impurity in Ranitidine Using the Agilent 6470 Triple Quadrupole LC/MS (5994-1668EN)
- 7. Determination of Nitrosamine Impurities Using the High-Resolution Agilent 6546 LC/Q-TOF (<u>5994-1372EN</u>)
- Simultaneous Determination of Eight Nitrosamine Impurities in Metformin Extended-Release Tablets Using the Agilent 6470 Triple Quadrupole LC/MS (5994-2533EN)

Agilent Chemical Standards

Looking for high quality standards for additional applications? <u>Search our</u> <u>catalog of 5000 chemical standards</u> to find the best standard to use in your method, application or workflow. Can't find the right catalog standard? Our website also has an automated custom standard quoting tool to help you develop the standard that fits your needs.

Agilent LC Selector Tool

Are you looking for an LC column best suited for your applications? Use Agilent's <u>LC column selector tool</u> to select the optimal LC column.

Agilent CrossLab: Real insight, real outcomes

CrossLab goes beyond instrumentation to bring you services, consumables, and lab-wide resource management. So your lab can improve efficiency, optimize operations, increase instrument uptime, develop user skill, and more.



Learn more about Agilent CrossLab, and see examples of insight that leads to great outcomes, at

www.agilent.com/crosslab

Find a local Agilent customer center in your country: www.agilent.com/chem/contactus

Get answers to your technical questions and access resources in the Agilent Community: **community.agilent.com**

U.S. and Canada 1-800-227-9770 agilent_inquiries@agilent.com

Europe info_agilent@agilent.com

Asia Pacific inquiry_lsca@agilent.com

This information is subject to change without notice.

RA44239.5764236111

© Agilent Technologies, Inc. 2021 Published in the USA, March 10, 2021 5994-2977EN

